

Product Datasheet

Mouse CLU ELISA Kit (orb775175)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	18.75-1200 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q5SW19
Sensitivity	7.3 ng/mL
Expiration Date	Please enquire.

