

Product Datasheet

Mouse FGF21 ELISA Kit (orb775249)

Description This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...

Reactivity Mouse

Range 3.13-200 pg/mL

Tested Applications ELISA

Note For research use only

Sample Types serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [Q9JJN1](#)

Sensitivity 5.2 pg/mL

Expiration Date Please enquire.

