

Product Datasheet

Mouse GIP ELISA Kit (orb775446)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P48756
Sensitivity	25.2 pg/mL
Expiration Date	Please enquire.

