

## Product Datasheet

### Mouse CASP9 ELISA Kit (orb775566)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates, cell culture supernates and other biological fluids

**Uniprot ID** [Q8C3Q9](#)

**Sensitivity** 0.072 ng/mL

**Expiration Date** Please enquire.

