

Product Datasheet

Human TFF2 ELISA Kit (orb775634)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q03403
Sensitivity	24.7 pg/mL
Expiration Date	Please enquire.

