

## Product Datasheet

### Human dSIP ELISA Kit (orb775810)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro... |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Range</b>               | 23.44-1500 pg/mL                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates and other biological fluids                                           |
| <b>Sensitivity</b>         | 7.15 pg/mL                                                                                              |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |

