

## Product Datasheet

### Human SOCS3 ELISA Kit (orb776008)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 15.63-1000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** **O14543**

**Sensitivity** 6.1 pg/mL

**Expiration Date** Please enquire.

