

Product Datasheet

Human GDF2 ELISA Kit (orb776026)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	62.5-4000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9UK05
Sensitivity	24.8 pg/mL
Expiration Date	Please enquire.

