

## Product Datasheet

### Human PRKAb1 ELISA Kit (orb776110)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.32-20 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates and other biological fluids.

**Uniprot ID** [Q9Y478](#)

**Sensitivity** 0.124 ng/mL

**Expiration Date** Please enquire.

