

Product Datasheet

Rat TGFb1 ELISA Kit (orb776175)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids
Uniprot ID	P17246
Sensitivity	6.4 pg/mL
Expiration Date	Please enquire.

