

Product Datasheet

Human CA125 ELISA Kit (orb776186)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	1.57-100 U/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8WXI7
Sensitivity	0.58 U/mL
Expiration Date	Please enquire.

