

Product Datasheet

Mouse EDN1 ELISA Kit (orb776261)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	7.82-500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q61614
Sensitivity	2.43 pg/mL
Expiration Date	Please enquire.

