

Product Datasheet

Rabbit MCP1 ELISA Kit (orb776287)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rabbit
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P28292
Sensitivity	5.4 pg/mL
Expiration Date	Please enquire.

