

Product Datasheet

Rat PROCR ELISA Kit (orb776341)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q4V8I1
Sensitivity	0.55 ng/mL
Expiration Date	Please enquire.

