


Product Datasheet

Human CCK ELISA Kit (orb776498)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P06307
Sensitivity	4.93 pg/mL
Expiration Date	Please enquire.



Standard Curve of Human CCK ELISA Kit.

