

Product Datasheet

Rabbit MMP13 ELISA Kit (orb776649)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rabbit
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O62806
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

