

Product Datasheet

Mouse GC ELISA Kit (orb776651)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Mouse
Range	23.44-1500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P55095
Sensitivity	7.22 pg/mL
Expiration Date	Please enquire.

