

Product Datasheet

Mouse FGL2 ELISA Kit (orb776661)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P12804
Sensitivity	0.135 ng/mL
Expiration Date	Please enquire.

