


Product Datasheet

OGDH antibody (orb77667)

Description	Rabbit polyclonal antibody to OGDH
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 1-288aa of human OGDH protein
Target	OGDH
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Dilution Range	WB: 1:500-1:4,000
Expiration Date	12 months from date of receipt.



Western blot analysis of fetal thymus ly...



Immunohistochemical analysis of formalin...