

Product Datasheet

Human FLNC ELISA Kit (orb776921)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q14315
Sensitivity	0.126 ng/mL
Expiration Date	Please enquire.

