

## Product Datasheet

### Mouse PSGL1 ELISA Kit (orb776952)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Serum, plasma and other biological fluids.

**Uniprot ID** [Q62170](#)

**Sensitivity** 0.066 ng/mL

**Expiration Date** Please enquire.

