

Product Datasheet

Mouse ANGPTL2 ELISA Kit (orb777053)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9R045
Sensitivity	12.5 pg/mL
Expiration Date	Please enquire.

