

## Product Datasheet

### Human PRDX2 ELISA Kit (orb777123)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 31.25-2000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** serum, plasma, tissue homogenates and other biological fluids

**Uniprot ID** **P32119**

**Sensitivity** 14.2 pg/mL

**Expiration Date** Please enquire.

