

Product Datasheet

Human TNNI1 ELISA Kit (orb777276)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	46.88-3000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P19237
Sensitivity	17.59 pg/mL
Expiration Date	Please enquire.

