

## Product Datasheet

### Human CCDC80 ELISA Kit (orb777369)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** serum, plasma, tissue homogenates and other biological fluids

**Uniprot ID** [Q76M96](#)

**Sensitivity** 0.054 ng/mL

**Expiration Date** Please enquire.

