

Product Datasheet

Mouse bACE1 ELISA Kit (orb777389)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cerebrospinal fluid and other biological fluids
Uniprot ID	P56818
Sensitivity	14.1 pg/mL
Expiration Date	Please enquire.

