


Product Datasheet

Human PLA2R1 ELISA Kit (orb777422)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids
Uniprot ID	Q13018
Sensitivity	0.056 ng/mL
Expiration Date	Please enquire.



Standard Curve of Human PLA2R1 ELISA Kit...

