

## Product Datasheet

### Mouse ABCG2 ELISA Kit (orb777443)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Mouse

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates and other biological fluids.

**Uniprot ID** [Q7TMS5](#)

**Sensitivity** 0.055 ng/mL

**Expiration Date** Please enquire.

