



Product Datasheet

ADAP2 antibody (orb77749)

Description	Rabbit polyclonal antibody to ADAP2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 48-281aa of human ADAP2 protein
Target	ADAP2
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Dilution Range	WB: 1:500-1:2,000. IHC: 1:100-1:200. ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of fetal spleen an...