

Product Datasheet

Rat PDCN ELISA Kit (orb777550)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	urine, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q8K4G9
Sensitivity	0.056 ng/mL
Expiration Date	Please enquire.

