

Product Datasheet

Mouse GDF2 ELISA Kit (orb777570)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9WV56
Sensitivity	28 pg/mL
Expiration Date	Please enquire.

