

Product Datasheet

Mouse ACE2 ELISA Kit (orb777589)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	12.5-800 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q8R0I0
Sensitivity	5.5 pg/mL
Expiration Date	Please enquire.

