

## Product Datasheet

### Human CHRDL1 ELISA Kit (orb777634)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.16-10 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids

**Uniprot ID** [Q9BU40](#)

**Sensitivity** 0.058 ng/mL

**Expiration Date** Please enquire.

