

## Product Datasheet

### Human CYSLTR1 ELISA Kit (orb777807)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.32-20 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids.

**Uniprot ID** [Q9Y271](#)

**Sensitivity** 0.135 ng/mL

**Expiration Date** Please enquire.

