

## Product Datasheet

### Human DDAH2 ELISA Kit (orb777959)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 0.32-20 ng/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Tissue homogenates, cell lysates and other biological fluids.

**Uniprot ID** **O95865**

**Sensitivity** 0.116 ng/mL

**Expiration Date** Please enquire.

