

Product Datasheet

Human MMRN2 ELISA Kit (orb777994)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	Q9H8L6
Sensitivity	29 pg/mL
Expiration Date	Please enquire.

