

## Product Datasheet

### Animal bLg ELISA Kit (orb778124)

|                            |                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...    |
| <b>Reactivity</b>          | Bovine                                                                                                     |
| <b>Range</b>               | 0.16-10 ng/mL                                                                                              |
| <b>Tested Applications</b> | ELISA                                                                                                      |
| <b>Note</b>                | For research use only                                                                                      |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, milk, cell culture supernates and other biological fluids |
| <b>Sensitivity</b>         | 0.067 ng/mL                                                                                                |
| <b>Expiration Date</b>     | Please enquire.                                                                                            |

