

Product Datasheet

Human TNFSF14 ELISA Kit (orb778342)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Human

Range 15.63-1000 pg/mL

Tested Applications ELISA

Note For research use only

Sample Types serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [O43557](#)

Sensitivity 5.5 pg/mL

Expiration Date Please enquire.

