

Product Datasheet

Rabbit PPAR γ ELISA Kit (orb778474)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rabbit
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	O19052
Sensitivity	0.28 ng/mL
Expiration Date	Please enquire.

