

## Product Datasheet

### Mouse CART ELISA Kit [Out of stock] (orb778677)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Mouse                                                                                                   |
| <b>Range</b>               | 0.16-10 ng/mL                                                                                           |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Sample Types</b>        | Tissue homogenates, cell lysates and other biological fluids.                                           |
| <b>Uniprot ID</b>          | <b>P56388</b>                                                                                           |
| <b>Sensitivity</b>         | 0.055 ng/mL                                                                                             |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |

