


Product Datasheet

IL6 antibody (orb77872)

Description	Rabbit polyclonal antibody to IL6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 30-212aa of human IL6 protein
Target	IL6
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	24kDa
Dilution Range	WB: 1:500-1:3,200
Expiration Date	12 months from date of receipt.



Western blot analysis of fetal spleen