



Product Datasheet

SCFD1 antibody (orb77875)

Description	Rabbit polyclonal antibody to SCFD1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 267-519aa of human SCFD1 protein
Target	SCFD1
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	72kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 and 293T cell lysate using SCFD1 antibody.



IHC analysis of HepG2 and 293T cell lysate using SCFD1 antibody.