

## Product Datasheet

### Human TMPO ELISA Kit (orb778784)

**Description** The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

**Reactivity** Human

**Range** 7.82-500 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** Serum, plasma, tissue homogenates and other biological fluids

**Uniprot ID** **P42166**

**Sensitivity** 12.9 pg/mL

**Expiration Date** Please enquire.

