

## Product Datasheet

### Human ADM2 ELISA Kit (orb778874)

**Description** This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...

**Reactivity** Human

**Range** 312.5-20000 pg/mL

**Tested Applications** ELISA

**Note** For research use only

**Sample Types** serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

**Uniprot ID** [Q7Z4H4](#)

**Sensitivity** 91.5 pg/mL

**Expiration Date** Please enquire.

