

Product Datasheet

Human TAC1 ELISA Kit (orb778892)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Human
Range	156.25-10000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P20366
Sensitivity	52.8 pg/mL
Expiration Date	Please enquire.

