

Product Datasheet

Rabbit S100B ELISA Kit (orb778949)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rabbit
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q6YNR6
Sensitivity	6.3 pg/mL
Expiration Date	Please enquire.

