

Product Datasheet

Human S100A16 ELISA Kit (orb778990)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	Q96FQ6
Sensitivity	0.069 ng/mL
Expiration Date	Please enquire.

