

Product Datasheet

Rat a1AGP ELISA Kit (orb779105)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	6.25-400 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P02764
Sensitivity	2.88 ng/mL
Expiration Date	Please enquire.

