

Product Datasheet

Rat APOC1 ELISA Kit (orb779164)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	1.57-100 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P19939
Sensitivity	0.64 ng/mL
Expiration Date	Please enquire.

