

Product Datasheet

Animal TGFb1 ELISA Kit (orb779316)

Description	This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate pro...
Reactivity	Bovine
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids
Uniprot ID	P18341
Sensitivity	4.67 pg/mL
Expiration Date	Please enquire.

